Barcelona, July, 20, 2020

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I, Laura López, Regulatory Affairs Manager at BioSystems S.A., hereby certify, that attached documents are true and accurate copies of User Manual of the product and translation of this document into Russian language in accordance with the documents.

Also, I hereby certify that the attached documents are for purposes of product registration in Russia and in the countries of Eurasian Economic Union.

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
Laura López
Regulatory Affairs Manager
BIOSYSTEMS S.A.

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BIOSYSTEMS, S.A.
Costa Brava, 30,
08030 BARCELONA
Tel. +34 933 11 00 00

Regulatory Affairs Department
BIOSYSTEMS S.A.
llopez@biosystems.es
Barcelona, July 20th, 2020

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B. Re
mi
W
to
dc

For f

STOP
Store
Comp
they
their
On br
the ar
Indica
"Test

ADDI
PROV
cod. '
C.

Bio
Hu

REAC
Reag
Worki
C (co
and fo

SAMP
Serun
Alanir
hepar

CALIF
A rea
month
chang

ALANINE AMINOTRANSFERASE (ALT/GPT)

COD 21533 8 x 60 mL + 8 x 15 mL	COD 23533 4 x 60 mL + 4 x 15 mL
Only for <i>in vitro</i> use in the clinical laboratory	

ALANINE AMINOTRANSFERASE (ALT/GPT) IFCC

INTENDED USE

Reagent for the measurement of alanine aminotransferase (ALT or GPT) concentration in human serum or plasma. The obtained values are useful as an aid in the diagnosis and monitoring of disorders of liver diseases, especially acute ones.

This reagent is for use in the BioSystems BA analyzers or in other analyzer with similar performance characteristics.

CLINICAL SIGNIFICANCE


The aminotransferases catalyze the formation of glutamic acid from 2-oxoglutarate by transfer of amino groups. ALT is normally present in various tissues but its higher concentrations are found in liver and kidney.

The serum concentration of ALT is elevated in hepatitis and other forms of hepatic disease associated with necrosis: infectious mononucleosis, cholestasis, cirrhosis, metastatic carcinoma of the liver, delirium tremens, and after administration of various drugs, such as opiates, salicylates or ampicillin^{1,2}. Serum ALT concentration can also be elevated in skeletal or cardiac muscle disease^{1,2}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE OF THE METHOD

Alanine aminotransferase (ALT or GPT) catalyzes the transfer of the amino group from alanine to 2-oxoglutarate, forming pyruvate and glutamate. The catalytic concentration is determined from the rate of decrease of NADH, measured at 340 nm, by means of the lactate dehydrogenase (LDH) coupled reaction³⁻⁶.

CONTENTS

	COD 21533	COD 23533
A. Reagent	8 x 60 mL	4 x 60 mL
B. Reagent	8 x 15 mL	4 x 15 mL

COMPOSITION

A. Reagent: Tris 150 mmol/L, L-alanine 750 mmol/L, lactate dehydrogenase > 1350 U/L, pH 7.3.

B. Reagent: NADH 1.9 mmol/L, 2-oxoglutarate 75 mmol/L, sodium hydroxide 148 mmol/L, sodium azide 9.5 g/L.

WARNING: H302: Harmful if swallowed. EUH031: Contact with acids liberates toxic gas. P301+P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. P330: Rinse mouth.

For further warnings and precautions, see the product safety data sheet (SDS).

STORAGE AND STABILITY

Store at 2-8 °C.

Components are stable once opened until the expiry date marked in the label if they are stored well closed and care is taken to prevent contamination during their use.

On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 2 months.

Indications of deterioration: Absorbance of the blank below the limit indicated in "Test Parameters".

ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

□ Pyridoxal phosphate (BioSystems cod. 11666):

C. Reagent: Pyridoxal phosphate 10 mmol/L. 5 mL.

□ Biochemistry Calibrator (BioSystems cod. 18011) or Biochemistry Calibrator Human (BioSystems cod. 18044).

REAGENT PREPARATION

Reagents are provided ready to use.

Working Reagent A with Pyridoxal Phosphate (Note 1): Pour 0.6 mL of Reagent C (cod 11666), into the Reagent A bottle. Mix gently. Stable for 6 days at 2-8°C and for 6 days in the refrigerated compartment of the analyzer.

SAMPLES

Serum and plasma collected by standard procedures.

Alanine aminotransferase in serum and plasma is stable for 7 days at 2-8°C. Use heparin or EDTA as anticoagulant⁸.

CALIBRATION

A reagent blank should be done every day and a calibration at least every 2 months (6 days for Working Reagent with Pyridoxal Phosphate), after reagent lot change or as required by quality control procedures.

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, 18009 and 18042) and II (cod. 18007, 18010 and 18043) to verify the accuracy of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits.

REFERENCE VALUES

Reaction temperature	37°C	30°C
Without pyr-P, up to ^{3,6}	41 U/L = 0.68 Ukat/L	29 U/L = 0.48 Ukat/L
With pyr-P, up to ^{3,4}	65 U/L = 1.08 Ukat/L	35 U/L = 0.58 Ukat/L

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

The metrological characteristics described below have been obtained using a BA400 analyzer and following the guidelines of the Clinical & Laboratory Standards Institute (CLSI).

□ Detection limit: 8.5 U/L = 0.14 Ukat/L.

□ Linearity limit: 500 U/L = 8.33 Ukat/L.

□ Precision:

Mean concentration	Repeatability (CV)	Within-laboratory (CV)
40.2 U/L = 0.67 Ukat/L	3.9 %	5.0 %
133 U/L = 2.21 Ukat/L	1.2 %	1.4 %

□ Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

LIMITATIONS OF THE PROCEDURE

□ Interferences: bilirubin (up to 20 mg/dL), hemolysis (hemoglobin up to 1000 mg/dL) and lipemia (triglycerides up to 200 mg/dL) do not interfere. Other drugs and substances may interfere⁷.

NOTES

1. The IFCC recommended method specifies the addition of pyridoxal phosphate.

BIBLIOGRAPHY

- Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACCC Press, 2001.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. Burtis CA, Ashwood ER, Brunis DE. WB Saunders Co, 2005.
- Sociedad Española de Química Clínica, Comité Científico, Comisión de Enzimas. Método recomendado para la determinación en rutina de la concentración catalítica de la alanina aminotransferasa en suero sanguíneo humano. *Quim Clin* 1987; 6: 241-244.
- IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 4. Reference procedure for the measurement of catalytic concentration of alanine aminotransferase. *Clin Chem Lab Med* 2002; 40: 718-724.
- IFCC reference procedures for measurement of catalytic concentrations of enzymes: corrigendum, notes and useful advice. *Clin Chem Lab Med* 2010; 48: 615-621.
- Gella FJ, Olivella T, Cruz Pastor M, Arenas J, Moreno R, Durban R and Gómez JA. A simple procedure for routine determination of aspartate aminotransferase and alanine aminotransferase with pyridoxal phosphate. *Clin Chim Acta* 1985; 153: 241-247.
- Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACCC Press, 2000.
- World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1, Rev.2; 2002.

TEST PARAMETERS

These reagents may be used in several automatic analyzers. Specific instructions for application in many of them are available on request.

R1: use Reagent A, or Reagent A with Pyridoxal Phosphate, R2: use Reagent B.

	BA200	BA400
GENERAL		
Name	ALT-GPT	ALT-GPT
Short name	ALT	ALT
Sample type	serum / plasma	serum / plasma
Analysis mode	kinetic bireagent	kinetic bireagent
Unit	U/L	U/L
Decimals	0	0
Reaction type	decreasing	decreasing
PROCEDURE		
Reading mode	monochromatic	monochromatic
Main filter	340	340
Reference filter	-	-
Sample	16.5	16.5
Vol. R1	160	160
Vol. R2	40	40
Reading 1 (cycle)	23	46
Reading 2 (cycle)	33	66
Predilution factor	-	-
CALIBRATION AND BLANK		
Blank type	distilled water	distilled water

Calibration mode	experimental calibrator	experimental calibrator	Kinetic blank limit	-	-
Number of calibrators	1	1	Linearity limit	500	500
Calibration curve	-	-	Substrate depletion	0.5	0.5

OPTIONS					
Blank absorbance limit	1.400	1.400			
M215331-05	BioSystems S.A. Costa Brava, 30.	08030 Barcelona (Spain)	03/2018	www.biosystems.es	

Quality System certified according to
EN ISO 13485 and EN ISO 9001 standards

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- Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. Burtis CA, Ashwood ER, Bruns DE. WB Saunders Co, 2005.
- Sociedad Española de Química Clínica, Comité Científico, Comisión de Enzimas. Método recomendado para la determinación en rutina de la concentración catalítica de la alanina aminotransferasa en suero sanguíneo humano. *Quim Clin* 1987; 6: 241-244.
- IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 4. Reference procedure for the measurement of catalytic concentration of alanine aminotransferase. *Clin Chem Lab Med* 2002; 40: 718-724.
- IFCC reference procedures for measurement of catalytic concentrations of enzymes; corrigendum, notes and useful advice. *Clin Chem Lab Med* 2010; 48: 615-621.
- Gella FJ, Olivella T, Cruz Pastor M, Arenas J, Moreno R, Durban R and Gómez JA. A simple procedure for routine determination of aspartate aminotransferase and alanine aminotransferase with pyridoxal phosphate. *Clin Chim Acta* 1985; 153: 241-247.
- Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
- World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1, Rev.2; 2002.

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