12/2019

Barcelona, July, 20, 2020

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FE3575198

I, Laura López, Regulatory Affairs Manager at BioSystems S.A., hereby certify, that attached documents are true and accurate copies of User Manual of the product and translation of this document into Russian language in accordance with the documents.

Also, I hereby certify that the attached documents are for purposes of product registration in Russia and in the countries of Eurasian Economic Union.

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Laura López  
Regulatory Affairs Manager  
BIOSYSTEMS S.A.

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COD

INTENDE

Reagent f  
obtained v  
evaluate th  
This reage

CLINICAL

Urea is sy  
in the urin  
Elevated  
protein ca  
or treatme  
Post-rena  
prostatic  
variability  
Clinical di  
both clinic

PRINCIPL

Urea in the  
can be me

CONTENTS

A. Reag  
B. Reag

COMPOSIT

A. Reagent:  
dehydroge  
B. Reagent: M  
WARNING:  
P301+P311:  
unwell. P33

STORAGE AM

Store at 2-8°C.  
Components ar  
well closed and  
On board stabili  
stable 30 days.  
Indications of  
Parameters".

WARNING AND

Exercise the nor  
available for prof  
with local guidel  
reported to BioS

ADDITIONAL M

Biochemistry Cal  
cod. 18044).

REAGENT PREF

Reagents are pro

SAMPLES

Serum, plasma or  
Urea in serum  
anticoagulant<sup>2</sup>  
Urea in urine is st

CALIBRATION

A reagent blank s  
lot change or as r

QUALITY CONTR

It is recommende  
cod. 18042) and II  
18054 and cod. 18  
Each laboratory s  
corrective action if

**UREA/BUN-UV**

COD 21516 8 x 60 mL + 8 x 15 mL	COD 23516 4 x 60 mL + 4 x 15 mL
Only for <i>in vitro</i> use in the clinical laboratory	

**BioSystems**

**UREA/BUN-UV**

UREASE/GLUTAMATE DEHYDROGENASE

**INTENDED USE**

Reagent for the measurement of urea concentration in human serum, plasma or urine. The obtained values are useful as an aid in the diagnosis and treatment of chronic renal failure, and to evaluate the function of renal glomerulus.

This reagent is for use in the BioSystems BA analyzers

**CLINICAL SIGNIFICANCE**

Urea is synthesized in the liver as a by-product of the deamination of amino acids. Its elimination in the urine represents the major route for nitrogen excretion.


Elevated urea concentration in plasma is found as a result of a high-protein diet, increased protein catabolism, after a gastrointestinal hemorrhage, mild dehydration, shock and heart failure or treatment with glucocorticoids (pre-renal uremia)<sup>1,2</sup>.

Post-renal uremia is caused by conditions that obstruct urine outflow: nephrolithiasis, tumor or prostatic hypertrophy. The usefulness of urea as an indicator of renal function is limited by the variability of its plasma concentration as a result of nonrenal factors<sup>1,2</sup>.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

**PRINCIPLE OF THE METHOD**

Urea in the sample consumes, by means of the coupled reactions described below, NADH that can be measured by spectrophotometry<sup>3,4</sup>.

**CONTENTS**

	COD 21516	COD 23516
A. Reagent	8 x 60 mL	4 x 60 mL
B. Reagent	8 x 15 mL	4 x 15 mL

**COMPOSITION**

A. Reagent: Tris 100 mmol/L, 2-oxoglutarate 5.6 mmol/L, urease > 140 U/mL, glutamate dehydrogenase > 140 U/mL, ethyleneglicol 220 g/L, sodium azide 0.95, pH 8.0.

B. Reagent: NADH 1.5 mmol/L, sodium azide 9.5 g/L

**WARNING: H302 Harmful if swallowed. EUH031: Contact with acids liberates toxic gas. P301+P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. P330: Rinse mouth.**

**STORAGE AND STABILITY**

Store at 2-8°C.

Components are stable once opened until the expiry date, marked in the label if they are stored well closed and care is taken to prevent contamination during their use.

On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 30 days.

Indications of deterioration: Absorbance of the blank below the limit indicated in \*Test Parameters\*.

**WARNING AND PRECAUTIONS**

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

**ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)**

Biochemistry Calibrator (BioSystems cod. 18011) or Biochemistry Calibrator Human (BioSystems cod. 18044)

**REAGENT PREPARATION**

Reagents are provided ready to use.

**SAMPLES**

Serum, plasma or urine collected by standard procedures.

Urea in serum or plasma is stable for 7 days at 2-8°C. Heparin is recommended as anticoagulant<sup>5</sup>.

Urea in urine is stable for 2 days at room temperature if microbial growth is prevented<sup>5</sup>.

**CALIBRATION**

A reagent blank should be done every day and a calibration at least every 30 days, after reagent lot change or as required by quality control procedures.

**QUALITY CONTROL**

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, cod. 18009 and cod. 18042) and II (cod. 18007, cod. 18010 and 18043) and the Biochemistry Control Urine (cod. 18054 and cod. 18066) to verify the accuracy of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits.

**REFERENCE VALUES**

Serum and plasma<sup>1</sup>: 12.8 - 42.8 mg/dL urea = 6 - 20 mg/dL BUN = 2.14 - 7.14 mmol/L urea. Concentrations in the neonatal period are lower, and in adults over 60 years of age are higher than in adults. Concentrations also tend to be slightly higher in males than in females.

Urine<sup>1</sup>: 26 - 43 g/24-h urea = 12 - 20 g/24 h BUN = 428 - 714 mmol/24-h urea

These ranges are given for orientation only, each laboratory should establish its own reference ranges.

**METROLOGICAL CHARACTERISTICS**

The metrological characteristics described below have been obtained using a BA400 analyzer and following the guidelines of the Clinical & Laboratory Standards Institute (CLSI).

- Detection limit: 3.69 mg/dL urea = 1.72 mg/dL BUN = 0.614 mmol/L urea

- Linearity limit: 300 mg/dL urea = 140 mg/dL BUN = 60 mmol/L urea. For samples with higher values, dilute manually or refer to the Test Parameterization for Automatic dilution (note that all these samples will be diluted with the same dilution ratio)

- Precision:

Serum Mean urea concentration	Repeatability (CV)	Within-laboratory (CV)
26.8 mg/dL = 4.47 mmol/L	3.5 %	5.0 %
137 mg/dL = 22.9 mmol/L	1.1 %	1.7 %
Urine Mean urea concentration	Repeatability (CV)	Within-laboratory (CV)
1291 mg/dL = 215 mmol/L	3.1 %	4.3 %
1771 mg/dL = 295 mmol/L	2.9 %	3.1 %

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

**LIMITATIONS OF THE PROCEDURE**

- Interferences: bilirubin (up to 30 mg/dL), hemolysis (hemoglobin up to 500 mg/dL) and lipemia (triglycerides up to 1625 mg/dL) do not interfere. Other drugs and substances may interfere<sup>6</sup>.

**BIBLIOGRAPHY**

1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Burtis CA, Ashwood ER, Bruns DE. WB Saunders Co, 2012
2. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACCPress, 2001
3. Talke H and Schubert GE. Enzymatische hamstoffbestimmung in blut und sern im optischen test nach Warbur. *Klinische Wochenschrift* 1965; 43: 174-175.
4. Gutmann I, Bergmeyer HU. Methods of enzymatic Analysis, ed Bergmeyer HU. Academic Press, NY, 1974, 4:1794-1798.
5. World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1, Rev 2, 2002
6. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACCPress, 2000.

**TEST PARAMETERS**

R1: use Reagent A

R2: use Reagent B

	BA200	BA400
<b>GENERAL</b>		
Name	UREA-BUN-UV	UREA-BUN-UV
Short name	UREA	UREA
Sample type	serum / plasma / urine	serum / plasma / urine
Analysis mode	fixed-time bireagent	fixed-time bireagent
Unit	mg/dL	mg/dL
Decimals	0	0
Reaction type	decreasing	decreasing
<b>PROCEDURE</b>		
Reading mode	monochromatic	monochromatic
Main filter	340	340
Reference filter	-	-
Sample	2	2
Vol. R1	240	240
Vol. R2	60	60
Reading 1 (cycle)	20	40
Reading 2 (cycle)	23	45
Predilution factor	1 / 50	1 / 50
Reduced factor	10	10
<b>CALIBRATION AND BLANK</b>		
Blank type	distilled water	distilled water
Calibration mode	experimental calibrator	experimental calibrator
Number of calibrators	1	1
Calibration curve	-	-
<b>OPTIONS</b>		
Blank absorbance limit	1.100	1.100
Kinetic blank limit	-	-
Linearity limit	300 / 300 / 15000	300 / 300 / 15000
Substrate depletion	-	-

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**DISPOSAL AND DESTRUCTION OF MEDICAL DEVICES**

National standards must be followed when transferring and disposing of the product, its components and packaging. Chemicals and finished and used reagents must be disposed of in accordance with regional regulations.

**LIABILITY AND GUARANTEES**

Any change to the set and/or exchange or replacement of any components of different batches from one set to another may negatively affect the expected results and the validity of the analysis as a whole. In the event of such changes and/or substitutions, any requirement to replace the set is invalid. Claims submitted due to incorrect interpretation of the analysis results by the client are also invalid.

However, in case of any claims, the manufacturer's liability should not exceed the cost of the kit. Damage caused to the kit during transportation is not the responsibility of the manufacturer. The manufacturer's warranties of BioSystems S. A. apply to the fact that this product meets the specifications stated on the label, but only if the product is used by qualified personnel in accordance with the instructions set out on the label and the materials are used before the expiration date indicated on the label.

**MAINTENANCE**

No requirement

**MAINTENANCE SERVICE**

Not required

**POTENTIAL CONSUMER**

Professional user: a laboratory specialist who has the appropriate work permit

**INFORMATION ABOUT THE PRESENCE OF A MEDICINAL PRODUCT IN A MEDICAL DEVICE**

Does not contain

**INDICATIONS FOR USE**

Not applicable

**CONTRAINDICATIONS TO USE**

Not applicable

**STERILITY**

It is delivered in non-sterile form. Sterility is not required

**INFORMATION ABOUT THE POSSIBILITY OF REUSE**

Not re-usable

**SYMBOLS and LABELS**

Symbol	Name
	See manual
	Catalogue number
	Medical device <i>in vitro</i>
	Lot
	European Community symbol
	Temperature
	Expire before
	Manufacturer

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**AUTORISED REPRESENTATIVE of the MANUFACTURER**

"MD-Consulting and Development" LLC, 129337, Russian Federation, Moscow, Yaroslavskoye Shosse, 6/2, office 106

**MANUFACTURER AND ADDRESS OF PRODUCTION**

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