04/2020
Barcelona, November, 20, 2020

FL0587850

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I, Laura López, Regulatory Affairs Manager at BioSystems S.A., hereby certify, that attached documents are true and accurate copies of User Manual of the product and translation of this document into Russian language in accordance with the documents.

Also, I hereby certify that the attached documents are for purposes of product registration in Russia and in the countries of Eurasian Economic Union.

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Laura López
Regulatory Affairs Manager
BIOSYSTEMS S.A.

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INTENDED USE
Reagent for the assessment of its v

CLINICAL BENEFIT
Serum lipase conc in amylase and lip period of time. Ele pancreatic duct by treatments with opi Based on clinical technologies and variations.

Clinical diagnosis both clinical and la

PRINCIPLE OF T
Lipase catalyzes glutaric acid-(6' intermediate, glu alkaline solution determined from

1,2-O-dilauryl-rac

glutaric ac

CONTENTS AN

- A. Reagent: 1 x calcium chlor
- B. Reagent: 1 x methylresor.

WARNING: H2
P210: Keep av vapours. P280: P403+P233 Stc

For further wan
S. Lipase Star on the label

Components fr antibodies anti cautiously as p

STORAGE AN

Store at 2-8°C
Components a well closed an
Indications of 0.80 at 580 nr

WARNING AI
Exercise the available for accordance v shall be repor

- ADDITIONAL**
- Biochemi (BioSyste
 - Analyzer, to read at

REAGENT R
Reagents ar

SAMPLES
Serum or pl
Lipase conc months at -

PROCEDU

1. Bring th
2. Pipette

Rea
Sam

3. Mix anc

COD 11760 1 x 30 mL

Only for *in vitro* use in the clinical laboratoryLIPASE
DGGR**INTENDED USE**

Reagent for the measurement of lipase concentration in human serum or plasma for the assessment of its variations in general population.

CLINICAL BENEFIT


Serum lipase concentration increases after an attack of acute pancreatitis. In general, increases in amylase and lipase run in parallel course, but the elevation of lipase persists for a longer period of time. Elevations in serum lipase concentration may be also due to obstruction of the pancreatic duct by a calculus or by carcinoma, in acute and chronic renal disease as well as in treatments with opiates^{1,2,3}.

Based on clinical guidelines and textbooks, and when used in conjunction with other diagnostic technologies and options, this medical information is useful for the assessment of Lipase variations.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE OF THE METHOD

Lipase catalyzes the hydrolysis of the chromogenic substrate 1,2-O-dilauryl-rac-glycerol-3-glutaric acid-(6'-methylresorufin)-ester to 1,2-O-dilauryl-rac-glycerol and an unstable intermediate, glutaric acid-(6'-methylresorufin)-ester. This decomposes spontaneously in alkaline solution to form glutaric acid and methylresorufin. The catalytic concentration is determined from the rate of the red dye formation measured at 560 nm^{1,4}.

**CONTENTS AND COMPOSITION**

A. Reagent: 1 x 20 mL. Bicine buffer 50 mmol/L, colipase \geq 1 mg/L, deoxycholate 1.6 mmol/L, calcium chloride 10 mmol/L, pH 8.0.

B. Reagent: 1 x 10 mL. Tartrate buffer 10 mmol/L, 1,2-O-dilauryl-rac-glycerol-3-glutaric acid-(6'-methylresorufin)-ester \geq 0.3 mmol/L, taurodesoxycholate 8.0 mmol/L, pH 4.0.

WARNING: H226: Flammable liquid and vapour. H317: May cause an allergic skin reaction. P210: Keep away from heat, sparks, open flames or hot surfaces. P261: Avoid breathing vapours. P280: Wear protective gloves, protective clothing, eye protection, face protection. P403+P233 Store in a well-ventilated place. Keep container tightly closed.

For further warnings and precautions, see the product safety data sheet (SDS).

S. Lipase Standard: 1 for 5 mL. Human lipase in human serum matrix. Concentration is given on the label.

Components from human origin have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for Hbs antigen. However, they should be handled cautiously as potentially infectious.

STORAGE AND STABILITY

Store at 2-8°C.

Components are stable once opened until the expiry date marked in the label if they are stored well closed and care is taken to prevent contamination during their use.

Indications of deterioration. Presence of particulate material or absorbance of the blank over 0.80 at 580 nm.

WARNING AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to Biosystems S.A.

ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

- Biochemistry Calibrator (BioSystems cod. 18011) or Biochemistry Calibrator Human (BioSystems cod. 18044).
- Analyzer, spectrophotometer or photometer with cell holder thermostatable at 37°C and able to read at 580 nm. Cuvettes with 1 cm light path.

REAGENT PREPARATION

Reagents are provided ready to use.

SAMPLES

Serum or plasma collected by standard procedures. Heparin may be used as anticoagulant. Lipase concentration in the sample is stable for 7 days at 20-25 °C, 21 days at 4-8°C and 12 months at -20°C⁵.

PROCEDURE

1. Bring the reagents and the instrument to reaction temperature.
2. Pipette into a cuvette (Note 1):

Reagent A	850 μ L
Sample/Standard	15 μ L

3. Mix and insert the cuvette into the instrument. Start the stopwatch.

4. After 1 minute, pipette into a cuvette:

Reagent B	500 μ L
-----------	-------------

5. Mix and after 3 minutes, record initial absorbance and at 1 minute intervals thereafter for 3 minutes.
6. Calculate the difference between consecutive absorbances, and the average absorbance difference per minute ($\Delta A/\text{min}$).

CALCULATIONS

The lipase concentration in the sample is calculated using the following general formula:

$$\frac{\Delta A/\text{min Sample}}{\Delta A/\text{min Standard}} \times C \text{ Standard} = U/L$$

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, 18009 and 18042) and II (cod. 18007, 18010 and 18043) to verify the accuracy of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

REFERENCE VALUES²

Serum or plasma: 13-60 U/L = 0.22-1.00 μ kat/L.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

- Detection limit: 4.89 U/L = 0.08 μ kat/L. Quantification limit: 8.52 U/L = 0.14 μ kat/L.
- Linearity limit: 250 U/L = 4.17 μ kat/L. For higher values dilute sample 1/2 with distilled water and repeat measurement. Measuring range: (8.52 U/L = 0.14 μ kat/L) - (250 U/L = 4.17 μ kat/L).
- Precision:

Mean concentration	Repeatability (CV)	Within-laboratory (CV)
52.1 U/L = 0.86 μ kat/L	1.4 %	4.1 %
67.2 U/L = 1.11 μ kat/L	2.2 %	5.3 %
122 U/L = 2.03 μ kat/L	1.6 %	5.0 %

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Lipase values for human serum and heparin plasma samples obtained on BA400 analyzer (y) were compared with those determined on a Roche Cobas 8000 analyzer (x). Serum: Sample size (n)=88; Linear regression $y=5.99+0.993x$, $r=0.971$. The sample concentrations were between 15 and 240 U/L. Plasma: Sample size (n)=57; Linear regression $y=5.59+0.964x$, $r=0.996$. The sample concentrations were between 12 and 285 U/L.

LIMITATIONS OF THE PROCEDURE

- Interferences: hemolysis (hemoglobin up to 500 mg/dL), bilirubin (up to 30 mg/dL) and lipemia (triglycerides up to 300 mg/dL) do not interfere. Other drugs and substances may interfere⁶.

The triglycerides reagent contains a very high lipase concentration that could interfere in lipase measurements. Current versions of Software solve the contamination by means of protocols already programmed by default in the analyzer. It is recommended to perform lipase measurements in series without triglycerides assays when abnormal results are detected.

NOTES

1. This reagent may be used in several automatic analyzers. Instructions for many of them are available upon request.

BIBLIOGRAPHY

1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 6th ed. Rifai N, Horvath AR, Wittwer CT. WB Saunders Co, 2018.
2. Junge W, Abicht K, Goldman J et al. Evaluation of the colorimetric liquid assay for pancreatic lipase on Hitachi analyzers in clinical centers in Europe, Japan, and USA. *Clin Chem Lab Med* 1999;37, special suppl:469.
3. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
4. Panteghini M, Bonora R, Pagani F. Measurement of pancreatic lipase activity in serum by a kinetic colorimetric assay using a new chromogenic substrate. *Ann Clin Biochem* 2001;38:365-370.
5. World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1, Rev.2; 2002.
6. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

COD 21760 1 x 20 mL + 1 x 10 mL
Only for <i>in vitro</i> use in the clinical laboratory

**INTENDED USE**

Reagent for the measurement of lipase concentration in human serum or plasma for the assessment of its variations in general population.

This reagent is for use in the BioSystems BA analyzers.

CLINICAL BENEFIT


Serum lipase concentration increases after an attack of acute pancreatitis. In general, increases in amylase and lipase run in parallel course, but the elevation of lipase persists for a longer period of time. Elevations in serum lipase concentration may be also due to obstruction of the pancreatic duct by a calculus or by carcinoma, in acute and chronic renal disease as well as in treatments with opiates^{1,2,3}.

Based on clinical guidelines and textbooks, and when used in conjunction with other diagnostic technologies and options, this medical information is useful for the assessment of Lipase variations.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE OF THE METHOD

Lipase catalyzes the hydrolysis of the chromogenic substrate 1,2-O-dilauryl-rac-glycerol-3-glutaric acid-(6'-methylresorufin)-ester to 1,2-O-dilauryl-rac-glycerol and an unstable intermediate, glutaric acid-(6'-methylresorufin)-ester. This descomposes spontaneously in alkaline solution to form glutaric acid and methylresorufin. The catalytic concentration is determined from the rate of the red dye formation measured at 560 nm^{1,4}.

**CONTENTS AND COMPOSITION**

A. Reagent: 1 x 20 mL. Bicine buffer 50 mmol/L, colipase \geq 1 mg/L, deoxycholate 1.6 mmol/L, calcium chloride 10 mmol/L, pH 8.0.

B. Reagent: 1 x 10 mL. Tartrate buffer 10 mmol/L, 1,2-O-dilauryl-rac-glycerol-3-glutaric acid-(6'-methylresorufin)-ester \geq 0.3 mmol/L, taurodesoxycholate 8.0 mmol/L, pH 4.0.

WARNING: H226: Flammable liquid and vapour. H317: May cause an allergic skin reaction. P210: Keep away from heat, sparks, open flames or hot surfaces. P261: Avoid breathing vapours. P280: Wear protective gloves, protective clothing, eye protection, face protection. P403+P233 Store in a well-ventilated place. Keep container tightly closed.

STORAGE AND STABILITY

Store at 2-8°C.

Components are stable once opened until the expiry date marked in the label if they are stored well closed and care is taken to prevent contamination during their use.

Indications of deterioration: Absorbance of the blank over the limit indicated in "Test Parameters".

On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 30 days.

WARNING AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to Biosystems S.A.

ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

Biochemistry Calibrator (BioSystems cod. 18011), or Biochemistry Calibrator Human (BioSystems cod. 18044).

REAGENT PREPARATION

Reagents are provided ready to use.

SAMPLES

Serum or plasma collected by standard procedures. Heparin may be used as anticoagulant.

Lipase concentration in the sample is stable for 7 days at 20-25 °C, 21 days at 4-8°C and 12 months at -20°C⁵.

CALIBRATION

A reagent blank should be done every day and a calibration at least every 30 days, after reagent lot change or as required by quality control procedures.

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, 18009 and 18042) and II (cod. 18007, 18010 and 18043) to verify the accuracy of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits.

REFERENCE VALUES²

Serum or plasma: 13-60 U/L = 0.22-1.00 μ kat/L.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

The metrological characteristics described below have been obtained using a BA400 analyzer and following the guidelines of the Clinical & Laboratory Standards Institute (CLSI).

- Detection limit: 4.89 U/L = 0.08 μ kat/L. Quantification limit: 9.52 U/L = 0.14 μ kat/L.
- Linearity limit: 250 U/L = 4.17 μ kat/L. For higher values dilute sample 1/2 with distilled water and repeat measurement. Measuring range: (8.52 U/L = 0.14 μ kat/L) - (250 U/L = 4.17 μ kat/L).
- Precision:

Mean concentration	Repeatability (CV)	Within-laboratory (CV)
52.1 U/L = 0.86 μ kat/L	1.4 %	4.1 %
67.2 U/L = 1.11 μ kat/L	2.2 %	5.3 %
122 U/L = 2.03 μ kat/L	1.6 %	5.0 %

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Lipase values for human serum and heparin plasma samples obtained on BA400 analyzer (y) were compared with those determined on a Roche Cobas 8000 analyzer (x). Serum: Sample size (n)=88; Linear regression $y=5.99+0.993x$, $r=0.971$. The sample concentrations were between 15 and 240 U/L. Plasma: Sample size (n)=57; Linear regression $y=5.59+0.964x$, $r=0.996$. The sample concentrations were between 12 and 285 U/L.

LIMITATIONS OF THE PROCEDURE

- Interferences: hemolysis (hemoglobin up to 500 mg/dL), bilirubin (up to 30 mg/dL) and lipemia (triglycerides up to 300 mg/dL) do not interfere. Other drugs and substances may interfere⁶.
- The triglycerides reagent contains a very high lipase concentration that interferes in lipase measurements by contamination of the reaction cuvette that is not eliminated with ordinary washing. It is recommended to perform lipase measurements in series without triglycerides assays to avoid reaction cuvette well contamination.

BIBLIOGRAPHY

1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 6th ed. Rifai N, Horvath AR, Wittwer CT. WB Saunders Co, 2018.
2. Junge W, Abicht K, Goldman J et al. Evaluation of the colorimetric liquid assay for pancreatic lipase on Hitachi analyzers in clinical centers in Europe, Japan, and USA. *Clin Chem Lab Med* 1999;37, special suppl:469.
3. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
4. Panteghini M, Bonora R, Pagani F. Measurement of pancreatic lipase activity in serum by a kinetic colorimetric assay using a new chromogenic substrate. *Ann Clin Biochem* 2001;38:365-370.
5. World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1, Rev.2; 2002.
6. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

TEST PARAMETERS

R1: use Reagent A

R2: use Reagent B

	BA200	BA400
GENERAL		
Name	LIP DGGR	LIP DGGR
Short name	LIP DGGR	LIP DGGR
Sample type	serum / plasma	serum / plasma
Analysis mode	kinetic bireagent	kinetic bireagent
Unit	U/L	U/L
Decimals	2	2
Reaction type	increasing	increasing
PROCEDURE		
Reading mode	monochromatic	monochromatic
Main filter	560	560
Reference filter	-	-
Sample	3	3
Vol. R1	170	170
Vol. R2	100	100
Reading 1 (cycle)	28	56
Reading 2 (cycle)	34	67
Predilution factor	-	-
Reduced factor	2	2
CALIBRATION AND BLANK		
Blank type	distilled water	distilled water
Calibration mode	experimental calibrator	experimental calibrator
Number of calibrators	1	1
Calibration curve	-	-
OPTIONS		
Blank absorbance limit	0.800	0.800
Kinetic blank limit	-	-
Linearity limit	250	250
Substrate depletion	-	-

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1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 6th ed. Rifai N, Horvath AR, Wittwer CT. WB Saunders Co, 2018.

2. Junge W, Abicht K, Goldman J et al. Evaluation of the colorimetric liquid assay for pancreatic lipase on Hitachi analyzers in clinical centers in Europe, Japan, and USA. Clin Chem Lab Med 1999;37, special suppl:469.

3. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.

4. Panteghini M, Bonora R, Pagani F. Measurement of pancreatic lipase activity in serum by a kinetic colorimetric assay using a new chromogenic substrate. Ann Clin Biochem 2001;38:365-370.

5. World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1, Rev.2; 2002.

6. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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