03/2022
Barcelona, November, 17, 2022

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I, Laura López, Regulatory Affairs Manager at BioSystems S.A., hereby certify, that attached documents are true and accurate copies of User Manual of the product and translation of this document into Russian language in accordance with the documents.

Also, I hereby certify that the attached documents are for purposes of product registration in Russia and in the countries of Eurasian Economic Union.

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Laura Lopez
Regulatory Affairs Manager
BIOSYSTEMS S.A.

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BIOSYSTEMS, S.A.
Costa Brava, 30,
08030 BARCELONA
Tel. +34 933 11 00 00

Regulatory Affairs Department
BIOSYSTEMS S.A.
llopez@biosystems.es
Barcelona, November, 17th, 2022

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COD 11534 1 x 40 mL
STORE AT 2-8°C
Reagents for measurement of α -amylase concentration Only for <i>in vitro</i> use in the clinical laboratory

α -AMYLASE-EPS

BioSystems

α -AMYLASE-EPS
IFCC

PRINCIPLE OF THE METHOD

Amylase catalyzes the hydrolysis of 4-nitrophenyl-maltoheptaoside-ethylidene to smaller oligosaccharides which are hydrolyzed by α -glucosidase liberating 4-nitrophenol. The catalytic concentration is determined from the rate of 4-nitrophenol formation, measured at 405 nm^{1,2,3}.

CONTENTS AND COMPOSITION

- A. Reagent. 1 x 32 mL. HEPES 50 mmol/L, calcium chloride 0.075 mmol/L, sodium chloride 90 mmol/L, magnesium chloride 13 mmol/L, α -glucosidase > 4 U/mL, pH 7.1.
B. Reagent. 1 x 8 mL. HEPES 50 mmol/L, 4-Nitrophenyl-maltoheptaoside-ethylidene 18 mmol/L, pH 7.1.

STORAGE

Store at 2-8°C.

Reagents are stable until the expiry date shown on the label when stored tightly closed and if contaminations are prevented during their use.

Indications of deterioration:

- Reagents: Presence of particulate material, turbidity, absorbance of the blank over 0.300 at 405 nm (1 cm cuvette).

WARNING AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

REAGENT PREPARATION

Working Reagent: Pour the contents of the Reagent B into the Reagent A bottle. Mix gently. Other volumes can be prepared in the proportion: 4 mL Reagent A + 1 mL Reagent B. Stable for 20 days at 2-8°C.

ADDITIONAL EQUIPMENT

- Analyzer, spectrophotometer or photometer with cell holder thermostatable at 37°C and able to read at 405 nm.
- Cuvettes with 1 cm light path.

SAMPLES

Serum, plasma or urine collected by standard procedures.

α -Amylase in serum or plasma is stable for 1 month at 2-8°C. Use heparin or EDTA as anticoagulant.

α -Amylase in urine is stable for 1 month at 2-8°C if pH is adjusted to approximately 7 before storage.

PROCEDURE

- Bring the Reagent and the instrument to reaction temperature.
- Pipette into a cuvette: (Notes 1,2)

	Serum or plasma 37°C	Urine 37°C
Working Reagent	1.0 mL	1.0 mL
Sample	30 μ L	15 μ L

- Mix and insert the cuvette into the photometer. Start the stopwatch.
- Record initial absorbance and at 1 minute intervals thereafter for 3 minutes.
- Calculate the difference between consecutive absorbances, and the average absorbance difference per minute ($\Delta A/min$).

CALCULATIONS

The α -Amylase concentration in the sample is calculated using the following general formula:

$$\Delta A/min \times \frac{Vt \times 10^6}{\epsilon \times l \times Vs} = U/L$$

The molar absorbance (ϵ) of 4-nitrophenol at 405 nm is 10600 and the lightpath (l) is 1 cm. For serum and plasma samples, the total reaction volume (Vt) is 1.030 at 37°C and the sample volume (Vs) is 0.030 at 37°C. For urine samples, the total reaction volume (Vt) is 1.015 at 37°C and the sample volume (Vs) is 0.015 at 37°C. 1 U/L are 0.0166 μ kat/L. The following formulas are deduced for the calculation of the catalytic concentration:

		37°C
$\Delta A/min$	Serum, plasma	x 3230 = U/L x 53.8 = μ kat/L
	Urine	x 8384 = U/L x 105.9 = μ kat/L

REFERENCE VALUES

Reaction temperature	Serum, plasma		Urine	
	U/L	μ kat/L	U/L	μ kat/L
37°C, up to ⁴	28 - 100	0.47 - 1.67	16 - 491	0.26 - 8.15

These ranges are given for orientation only, each laboratory should establish its own reference ranges.

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, cod. 18009 and cod. 18042), level II (cod. 18007, cod. 18010 and cod. 18043) and the Biochemistry Control Urine (cod. 18054 and cod. 18066) to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

METROLOGICAL CHARACTERISTICS

The metrological characteristics described below have been obtained using a A25 analyzer and following the guidelines of the Clinical & Laboratory Standards Institute (CLSI).

- Detection limit: 2.79 U/L = 0.046 μ kat/L. Quantification limit: 13.0 U/L = 0.217 μ kat/L.
- Linearity limit: 1300 U/L = 21.7 μ kat/L. Measuring range: 13.0 - 1300 U/L (serum) = 28.9 - 2600 U/L (urine).
- Precision:

Serum. Mean concentration	Repeatability (CV)	Within-laboratory (CV)
55.0 U/L = 0.91 μ kat/L	2.4 %	3.8 %
211 U/L = 3.50 μ kat/L	1.4 %	2.1 %
1049 U/L = 17.4 μ kat/L	1.2 %	1.4 %

Urine. Mean concentration	Repeatability (CV)	Within-laboratory (CV)
182 U/L = 3.03 μ kat/L	1.3 %	2.3 %
418 U/L = 6.94 μ kat/L	1.5 %	2.3 %
1982 U/L = 32.9 μ kat/L	1.0 %	1.7 %

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

LIMITATIONS OF THE PROCEDURE

- Interferences: bilirubin (up to 30 mg/dL), hemolysis (hemoglobin up to 604 mg/dL) and lipemia (triglycerides up to 659 mg/dL) do not interfere. Other drugs and substances may interfere⁵.

DIAGNOSTIC CHARACTERISTICS

α -Amylase catalyzes the hydrolysis of α -1,4-linkages of carbohydrates constituted of α -D-glucose units. The result is the formation of dextrans, maltose and some glucose molecules. α -Amylase is produced mainly by the exocrine pancreas (P-type) and the salivary glands (S-type) but it is also found in other tissues.

Assays of amylase activity in serum and urine are largely of use in the diagnosis of pancreatic diseases such as acute or chronic pancreatitis. Hyperamylasemia can also be due to renal insufficiency, acute pain of the abdomen, tumors of the lungs and the ovaries, salivary glands lesions, macroamylasemia, diabetic ketoacidosis, biliary tract disease, cerebral trauma, chronic alcoholism and drugs (opiates)^{6,7}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

NOTES

- Saliva and skin do contain α -amylase, therefore never pipette by mouth and avoid skin contact with the reagents.
- This reagent may be used in several automatic analysers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 8. Reference procedure for the measurement of catalytic concentration of α -amylase. *Clin Chem Lab Med* 2006; 44: 1146-1155.
- IFCC reference procedures for measurement of catalytic concentrations of enzymes: corrigendum, notes and useful advice. *Clin Chem Lab Med* 2010; 48: 615-621.
- Lorenz K. Routine α -amylase assay using protected 4-nitrophenyl-1,4- α -D-maltoheptaoside and a novel α -glucosidase. *Clin Chem* 2000; 46: 644-649.
- Junga W, Werner W, Wilke B et al. Development and evaluation of assays for the determination of total and pancreatic amylase at 37°C according to the principle recommended by the IFCC. *Clin Biochem* 2001; 34: 607-615.
- Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th edition. Burtis CA, Ashwood ER. WB Saunders Co., 2005.
- Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.

COD 12535 1 x 40 mL
Only for <i>in vitro</i> use in the clinical laboratory

INTENDED USE

Reagent for the measurement of α-amylase concentration in human serum, plasma or urine. The obtained values are useful as an aid in the diagnosis and treatment of acute and chronic pancreatitis.

These reagents are for use in the BioSystems A25 and A15 analyzer or in other analyzer with similar performance characteristics.

CLINICAL SIGNIFICANCE

α-Amylase catalyzes the hydrolysis of α-1,4-linkages of carbohydrates constituted of α-D-glucose units. The result is the formation of dextrins, maltose and some glucose molecules. α-Amylase is produced mainly by the exocrine pancreas (P-type) and the salivary glands (S-type) but it is also found in other tissues.

Assays of amylase activity in serum and urine are largely of use in the diagnosis of pancreatic diseases such as acute or chronic pancreatitis. Hyperamylasemia can also be due to renal insufficiency, acute pain of the abdomen, tumors of the lungs and the ovaries, salivary glands lesions, macroamylasemia, diabetic ketoacidosis, biliary tract disease, cerebral trauma, chronic alcoholism and drugs (opiates)^{1,2}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE OF THE METHOD

Amylase catalyzes the hydrolysis of 4-nitrophenyl-maltoheptaoside-ethylidene to smaller oligosaccharides which are hydrolyzed by α-glucosidase liberating 4-nitrophenol. The catalytic concentration is determined from the rate of 4-nitrophenol formation, measured at 405 nm^{3,4}.

CONTENTS AND COMPOSITION

- A. Reagent: 1 x 32 mL. HEPES 50 mmol/L, calcium chloride 0.075 mmol/L, sodium chloride 90 mmol/L, magnesium chloride 13 mmol/L, α-glucosidase > 4 U/mL, pH 7.1.
- B. Reagent: 1 x 8 mL. HEPES 50 mmol/L, 4-Nitrophenyl-maltoheptaoside-ethylidene 18 mmol/L, pH 7.1.

STORAGE AND STABILITY

Store at 2-8 °C.
Components are stable once opened until the expiry date marked in the label if they are stored well closed and care is taken to prevent contamination during their use.
On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 20 days.
Indications of deterioration: Absorbance of the blank over the limit indicated in "Test Parameters".

ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

Biochemistry Calibrator (BioSystems cod. 18011) or Biochemistry Calibrator Human (BioSystems cod. 18044).

REAGENT PREPARATION

Working Reagent: Pour the contents of the Reagent B into the Reagent A bottle. Mix gently. Other volumes can be prepared in the proportion: 4 mL Reagent A + 1 mL Reagent B. Stable for 20 days at 2-8°C.

SAMPLES

Serum, plasma or urine collected by standard procedures.
α-Amylase in serum or plasma is stable for 30 days at 2-8°C. Use heparin or EDTA as anticoagulant.
α-Amylase in urine is stable for 30 days at 2-8°C if pH is adjusted to approximately 7 before storage. Centrifuge or filter before testing.

CALIBRATION

A reagent blank should be done every day and a calibration at least every 20 days, after reagent lot change or as required by quality control procedures.

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, cod. 18009 and cod. 18042) and II (cod. 18007, cod. 18010 and cod. 18043) and the Biochemistry Control Urine (cod. 18054 and cod. 18066) to verify the accuracy of the measurement procedure.
Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits.

REFERENCE VALUES

Reaction temperature	Serum, plasma		Urine	
	U/L	μkat/L	U/L	μkat/L
30°C, up to ³	25-85	0.41-1.08	-	-
37°C, up to ⁶	28-100	0.47-1.67	18-491	0.28-8.15

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

The metrological characteristics described below have been obtained using a A25 analyzer and following the guidelines of the Clinical & Laboratory Standards Institute (CLSI).

- Detection limit: 2.79 U/L = 0.046 μkat/L. Quantification limit: 13.0 U/L = 0.217 μkat/L
- Linearity limit: 1300 U/L = 21.7 μkat/L. Measuring range: 13.0 - 1300 U/L (serum) = 28.9 - 2600 U/L (urine).
- Precision:

Serum. Mean concentration	Repeatability (CV)	Within-laboratory (CV)
55.0 U/L = 0.91 μkat/L	2.4 %	3.8 %
211 U/L = 3.50 μkat/L	1.4 %	2.1 %
1049 U/L = 17.4 μkat/L	1.2 %	1.4 %

Urine. Mean concentration	Repeatability (CV)	Within-laboratory (CV)
182 U/L = 3.03 μkat/L	1.3 %	2.3 %
418 U/L = 6.94 μkat/L	1.5 %	2.3 %
1982 U/L = 32.9 μkat/L	1.0 %	1.7 %

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

LIMITATIONS OF THE PROCEDURE

- Interferences: bilirubin (up to 30 mg/dL); hemolysis (hemoglobin up to 604 mg/dL) and lipemia (triglycerides up to 659 mg/dL) do not interfere. Other drugs and substances may interfere⁷.

BIBLIOGRAPHY

1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th edition. Burtis CA, Ashwood ER. WB Saunders Co., 2005.
2. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
3. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 8. Reference procedure for the measurement of catalytic concentration of α-amylase. *Clin Chem Lab Med* 2006; 44: 1146-1155.
4. IFCC reference procedures for measurement of catalytic concentrations of enzymes: corrigendum, notes and useful advice. *Clin Chem Lab Med* 2010, 48: 615-621.
5. Lorentz K. Routine α-amylase assay using protected 4-nitrophenyl-1,4-α-D-maltoheptaoside and a novel α-glucosidase. *Clin Chem* 2000;46:644-649.
6. Junge W, Werner W, Wilke B et al. Development and evaluation of assays for the determination of total and pancreatic amylase at 37°C according to the principle recommended by the IFCC. *Clin Biochem* 2001;34:607-615.
7. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

TEST PARAMETERS

	A25	A15
GENERAL		
Name	AMYLASE EPS	AMYLASE EPS
Short name	AMY EPS	AMY EPS
Sample type	SER / URI	SER / URI
Analysis mode	monoreagent kinetics	monoreagent kinetics
Unit	U/L	U/L
Decimals	0	0
Reaction type	increasing	increasing
PROCEDURE		
Reading mode	monochromatic	monochromatic
Main filter	405	405
Reference filter	-	-
Sample	9	9
Vol. R1	300	300
Vol. R2	-	-
Reading 1 (cycle)	5	4
Reading 2 (cycle)	14	10
Predilution factor	- / 2	- / 2
CALIBRATION AND BLANK		
Calibration mode	multiple	multiple
Number of calibrators	-	-
Calibration curve	-	-
OPTIONS		
Blank absorbance limit	0.300	0.300
Kinetic blank limit	-	-
Linearity limit	1300 / 2600	1300 / 2600
Substrate depletion	-	-

COD 21534 2 x 60 mL + 2 x 15 mL

Only for *in vitro* use in the clinical laboratory

INTENDED USE

Reagent for the measurement of α-amylase concentration in human serum, plasma or urine. The obtained values are useful as an aid in the diagnosis and treatment of acute and chronic pancreatitis.

This reagent is for use in the BioSystems BA analyzers or in other analyzer with similar performance characteristics.

CLINICAL SIGNIFICANCE

α-Amylase catalyzes the hydrolysis of α-1,4-linkages of carbohydrates constituted of α-D-glucose units. The result is the formation of dextrans, maltose and some glucose molecules. α-Amylase is produced mainly by the exocrine pancreas (P-type) and the salivary glands (S-type) but it is also found in other tissues.

Assays of amylase activity in serum and urine are largely of use in the diagnosis of pancreatic diseases such as acute or chronic pancreatitis. Hyperamylasemia can also be due to renal insufficiency, acute pain of the abdomen, tumors of the lungs and the ovaries, salivary glands lesions, macroamylasemia, diabetic ketoacidosis, biliary tract disease, cerebral trauma, chronic alcoholism and drugs (opiates)¹.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE OF THE METHOD

Amylase catalyzes the hydrolysis of 4-nitrophenyl-maltoheptaoside-ethylidene to smaller oligosaccharides which are hydrolyzed by α-glucosidase liberating 4-nitrophenol. The catalytic concentration is determined from the rate of 4-nitrophenol formation, measured at 405 nm³⁵.

CONTENTS AND COMPOSITION

A. Reagent: 2 x 60 mL. HEPES 50 mmol/L, calcium chloride 0.075 mmol/L, sodium chloride 90 mmol/L, magnesium chloride 13 mmol/L, α-glucosidase > 4 U/mL, pH 7.1.

B. Reagent: 2 x 15 mL. HEPES 50 mmol/L, 4-Nitrophenyl-maltoheptaoside-ethylidene 18 mmol/L, pH 7.1.

STORAGE AND STABILITY

Store at 2-8 °C.

Components are stable once opened until the expiry date marked in the label if they are stored well closed and care is taken to prevent contamination during their use.

On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 2 months.

Indications of deterioration: Absorbance of the blank over the limit indicated in "Test Parameters".

ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

Biochemistry Calibrator (BioSystems cod. 18011) or Biochemistry Calibrator Human (BioSystems cod. 18044).

REAGENT PREPARATION

Reagents are provided ready to use.

SAMPLES

Serum, plasma or urine collected by standard procedures.

α-Amylase in serum or plasma is stable for 30 days at 2-8°C. Use heparin or EDTA as anticoagulant.

α-Amylase in urine is stable for 30 days at 2-8°C if pH is adjusted to approximately 7 before storage. Centrifuge or filter before testing.

CALIBRATION

A reagent blank should be done every day and a calibration at least every 2 months, after reagent lot change or as required by quality control procedures.

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, cod. 18009 and cod. 18042) and II (cod. 18007, cod. 18010 and cod. 18043) and the Biochemistry Control Urine (cod. 18054 and cod. 18066) to verify the accuracy of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits.

REFERENCE VALUES

Reaction temperature	Serum, plasma		Urine	
	U/L	μkat/L	U/L	μkat/L
30°C, up to³	25-65	0.41-1.08	-	-
37°C, up to⁴	28-100	0.47-1.87	16-491	0.28-8.15

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

The metrological characteristics described below have been obtained using a BA400 analyzer and following the guidelines of the Clinical & Laboratory Standards Institute (CLSI).

- Detection limit: 5.6 U/L = 0.094 μkat/L.
- Linearity limit: 1300 U/L = 21.6 μkat/L.
- Precision:

Serum. Mean concentration	Repeatability (CV)	Within-laboratory (CV)
100 U/L = 1.67 μkat/L	1.5 %	1.9 %
203 U/L = 3.4 μkat/L	2.1 %	2.3 %

Urine. Mean concentration	Repeatability (CV)	Within-laboratory (CV)
103 U/L = 1.71 μkat/L	2.2 %	2.7 %
206 U/L = 3.42 μkat/L	2.8 %	3.1 %

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

LIMITATIONS OF THE PROCEDURE

- Interferences: bilirubin (up to 20 mg/dL), hemolysis (hemoglobin up to 1000 mg/dL) and lipemia (triglycerides up to 1000 mg/dL) do not interfere. Other drugs and substances may interfere⁷.

BIBLIOGRAPHY

1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th edition. Buris CA, Ashwood ER. WB Saunders Co., 2005.
2. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
3. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 8. Reference procedure for the measurement of catalytic concentration of α-amylase. *Clin Chem Lab Med* 2006; 44: 1146-1155.
4. IFCC reference procedures for measurement of catalytic concentrations of enzymes: compendium, notes and useful advice. *Clin Chem Lab Med* 2010; 48: 615-621.
5. Lorentz K. Routine α-amylase assay using protected 4-nitrophenyl-1,4-α-D-maltoheptaoside and a novel α-glucosidase. *Clin Chem* 2000;46:644-649.
6. Junge W, Werner W, Wilke B et al. Development and evaluation of assays for the determination of total and pancreatic amylase at 37°C according to the principle recommended by the IFCC. *Clin Biochem* 2001;34:607-615.
7. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

TEST PARAMETERS

These reagents may be used in several automatic analyzers. Specific instructions for application in many of them are available on request.

R1: use Reagent A, R2: use Reagent B.

	BA200	BA400
GENERAL		
Name	AMYLAASE EPS	AMYLAASE EPS
Short name	AMY EPS	AMY EPS
Sample type	serum / plasma / urine	serum / plasma / urine
Analysis mode	kinetic bireagent	kinetic bireagent
Unit	U/L	U/L
Decimals	0	0
Reaction type	increasing	increasing
PROCEDURE		
Reading mode	monochromatic	monochromatic
Main filter	405	405
Reference filter	-	-
Sample	6	9
Vol. R1	240	240
Vol. R2	60	60
Reading 1 (cycle)	21	43
Reading 2 (cycle)	31	63
Predilution factor	- / - / 2	- / - / 2
CALIBRATION AND BLANK		
Blank type	distilled water	distilled water
Calibration mode	experimental calibrator	experimental calibrator
Number of calibrators	1	1
Calibration curve	-	-
OPTIONS		
Blank absorbance limit	0,300	0,300
Kinetic blank limit	-	-
Linearity limit	1300 / 1300 / 2600	1300 / 1300 / 2600
Substrate depletion	-	-

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1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th edition Burtis CA, Ashwood ER. WB Saunders Co., 2005.
2. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AAC Press, 2001.
3. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C Part 8 Reference procedure for the measurement of catalytic concentration of α -amylase Clin Chem Lab Med 2006, 44: 1146-1155.
4. IFCC reference procedures for measurement of catalytic concentrations of enzymes: corrigendum, notes and useful advice. Clin Chem Lab Med 2010; 48: 615-621.
5. Lorentz K. Routine α -amylase assay using protected 4-nitrophenyl-1,4- β -D-maltoheptaoside and a novel α -glucosidase. Clin Chem 2000;46:644-649.
6. Junge W, Werner W, Wilke B et al. Development and evaluation of assays for the determination of total and pancreatic amylase at 37°C according to the principle recommended by the IFCC. Clin Biochem 2001;34:607615.
7. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AAC Press, 2000.

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1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th edition. Burtis CA, Ashwood ER. WB Saunders Co., 2005
2. Friedman and Young Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
3. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 8. Reference procedure for the measurement of catalytic concentration of α -amylase. Clin Chem Lab Med 2006; 44: 1146-1155.
4. IFCC reference procedures for measurement of catalytic concentrations of enzymes: corrigendum, notes and useful advice. Clin Chem Lab Med 2010; 48: 615-621.
5. Lorentz K. Routine α -amylase assay using protected 4-nitrophenyl-1,4- α -D-maltoheptaoside and a novel α -glucosidase. Clin Chem 2000;46:644-649.
6. Junge W, Werner W, Wilke B et al. Development and evaluation of assays for the determination of total and pancreatic amylase at 37°C according to the principle recommended by the IFCC. Clin Biochem 2001;34:607-615.
7. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

M21550r-07 BioSystems S.A. Costa Brava, 30. 08030 Barcelona (Spain)

www.biosystems.es

Quality System certified according to
EN ISO 13485 and EN ISO 9001 standards

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- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th edition Burtis CA, Ashwood ER. WB Saunders Co., 2005
- Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
- IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C Part 8 Reference procedure for the measurement of catalytic concentration of α -amylase. Clin Chem Lab Med 2006; 44: 1146-1155
- IFCC reference procedures for measurement of catalytic concentrations of enzymes: corrigendum, notes and useful advice. Clin Chem Lab Med 2010; 48: 615-621.
- Lorentz K. Routine α -amylase assay using protected 4-nitrophenyl-1,4- β -D-maltoheptaoside and a novel α -glucosidase. Clin Chem 2000; 46: 644-649
- Junge W, Werner W, Wilke B et al. Development and evaluation of assays for the determination of total and pancreatic amylase at 37°C according to the principle recommended by the IFCC. Clin Biochem 2001; 34: 607-615
- Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

M21550-07 BioSystems S.A. Costa Brava, 30. 08030 Barcelona (Spain)

www.biosystems.es

Quality System certified according to

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