12/2019

Barcelona, July, 20, 2020

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I, Laura López, Regulatory Affairs Manager at BioSystems S.A., hereby certify, that attached documents are true and accurate copies of User Manual of the product and translation of this document into Russian language in accordance with the documents.

Also, I hereby certify that the attached documents are for purposes of product registration in Russia and in the countries of Eurasian Economic Union.

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Laura López
Regulatory Affairs Manager
BIOSYSTEMS S.A.

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BIOSYSTEMS, S.A.
Costa Brava, 30,
08030 BARCELONA
Tel. +34 933 11 00 00

Regulatory Affairs Department
BIOSYSTEMS S.A.
llopez@biosystems.es
Barcelona, July 20th, 2020

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Clinical dia should inte

PRINCIPL

Aspartate group from catalytic c measured reaction³⁻⁶.

Oxalacetate

CONTENT

A. Reag

B. Reag

COMPOSIT

A. Reagent: 460 U/L, WARNING P280: V protection several m rinsing. P

3. Reagent: mmol/L, s WARNING toxic gas doctor/ph

or further w

TORAGE A

store at 2-8 :omponents ey are stor heir use.

In board sta re analyzer indications of Test Param

DDITIONAL ROVIDED)

od. 11666): C. Reager

Biochemist Human (Bi

EAGENT PI

agents are

ASPARTATE AMINOTRANSFERASE (AST/GOT)

COD 21531	8 x 60 mL + 8 x 15 mL	COD 23531	4 x 60 mL + 4 x 15 mL
Only for <i>in vitro</i> use in the clinical laboratory			

ASPARTATE AMINOTRANSFERASE (AST/GOT) IFCC

Working Reagent A with Pyridoxal Phosphate (Note 1): Pour 0.6 mL of Reagent C (cod 11666). Into the Reagent A bottle. Mix gently. Stable for 6 days at 2-8°C and for 6 days in the refrigerated compartment of the analyzer.

SAMPLES

Serum and plasma collected by standard procedures. Aspartate aminotransferase in serum and plasma is stable for 7 days at 2-8°C. Use heparin as anticoagulant⁸.

CALIBRATION

A reagent blank should be done every day and a calibration at least every 2 months (6 days for Working Reagent with Pyridoxal Phosphate), after reagent lot change or as required by quality control procedures.

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, 18009 and 18042) and II (cod. 18007, 18010 and 18043) to verify the accuracy of the measurement procedure. Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits.

REFERENCE VALUES

Reaction temperature	37°C	30°C
Without pyr-P, up to ¹	40 U/L = 0.67 Ukat/L	25 U/L = 0.42 Ukat/L
With pyr-P, up to ^{3,4}	50 U/L = 0.83 Ukat/L	30 U/L = 0.50 Ukat/L

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

The metrological characteristics described below have been obtained using a BA400 analyzer and following the guidelines of the Clinical & Laboratory Standards Institute (CLSI).

- Detection limit: 7.15 U/L = 0.119 Ukat/L.
- Linearity limit: 500 U/L = 8.33 Ukat/L.
- Precision:

Mean concentration	Repeatability (CV)	Within-laboratory (CV)
41.5 U/L = 0.69 Ukat/L	2.6 %	5.8 %
154 U/L = 2.55 Ukat/L	1.0 %	2.7 %

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

LIMITATIONS OF THE PROCEDURE

- Interferences: Lipemia (triglycerides 2 g/L) interfere. Bilirubin (20 mg/dL) and Hemolysis (hemoglobin 10 g/L) do not interfere. Other drugs and substances may interfere⁷.

NOTES

1. The IFCC recommended method specifies the addition of pyridoxal phosphate.

BIBLIOGRAPHY

1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. Burtis CA, Ashwood ER, Bruns DE. WB Saunders Co, 2005.
2. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AAC Press, 2001.
3. Sociedad Española de Química Clínica, Comité Científico, Comisión de Enzimas. Método recomendado para la determinación en rutina de la concentración catalítica de la aspartato aminotransferasa en suero sanguíneo humano. Quim Clin 1987; 6: 235-239.
4. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C, Part 5. Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 2002; 40:725-733.
5. IFCC reference procedures for measurement of catalytic concentrations of enzymes: corrigendum, notes and useful advice. Clin Chem Lab Med 2010; 48: 615-621.
6. Gella FJ, Olivella T, Cruz Pastor M, Areñas J, Moreno R, Durban R and Gómez JA. A simple procedure for routine determination of aspartate aminotransferase and alanine aminotransferase with pyridoxal phosphate. Clin Chim Acta 1985; 153: 241-247.
7. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AAC Press, 2000.
8. World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1, Rev.2, 2002.

TEST PARAMETERS

These reagents may be used in several automatic analyzers. Specific instructions for application in many of them are available on request.

R1: use Reagent A, or Reagent A with Pyridoxal Phosphate, R2: use Reagent B.

GENERAL	BA200	BA400
Name	AST-GOT	AST-GOT
Short name	AST/GOT	AST/GOT
Sample type	serum /plasma	serum /plasma
Analysis mode	kinetic bireagent	kinetic bireagent
Unit	U/L	U/L
Decimals	0	0
Reaction type	decreasing	decreasing

INTENDED USE

Reagent for the measurement of aspartate aminotransferase (AST or GOT) concentration in human serum or plasma. The obtained values are useful as an aid in the diagnosis and monitoring of disorders of liver diseases, especially acute ones.

This reagent is for use in the BioSystems BA analyzers or in other analyzer with similar performance characteristics.

CLINICAL SIGNIFICANCE

The aminotransferases catalyze the formation of glutamic acid from 2-oxoglutarate by transfer of amino groups. AST is found in highest concentration in the liver and heart muscle but it is also abundant in skeletal muscle, kidney and pancreas.


The serum concentration of AST is elevated in hepatitis and other forms of hepatic disease associated with necrosis: infectious mononucleosis, cholestasis, cirrhosis, metastatic carcinoma of the liver, delirium tremens, and after administration of various drugs^{1,2}.

Serum AST concentration is also elevated after myocardial infarction, in skeletal muscle disease (as progressive muscular dystrophy), in acute pancreatitis or hemolytic disease and other^{1,2}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE OF THE METHOD

Aspartate aminotransferase (AST or GOT) catalyzes the transfer of the amino group from aspartate to 2-oxoglutarate, forming oxalacetate and glutamate. The catalytic concentration is determined from the rate of decrease of NADH, measured at 340 nm, by means of the malate dehydrogenase (MDH) coupled reaction³⁻⁶.

CONTENTS

	COD 21531	COD 23531
A. Reagent	8 x 60 mL	4 x 60 mL
B. Reagent	8 x 15 mL	4 x 15 mL

COMPOSITION

A. Reagent: Tris 121 mmol/L, L-aspartate 362 mmol/L, malate dehydrogenase 460 U/L, lactate dehydrogenase 660 U/L, pH 7.8.

WARNING: H315: Causes skin irritation. H319: Causes serious eye irritation. P280: Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P332+P313: If skin irritation occurs: Get medical advice/attention.

B. Reagent: NADH 1.9 mmol/L, 2-oxoglutarate 75 mmol/L, sodium hydroxide 148 mmol/L, sodium azide 9.5 g/L.

WARNING: H302: Harmful if swallowed. EUH031: Contact with acids liberates toxic gas. P301+P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. P330: Rinse mouth.

For further warnings and precautions, see the product safety data sheet (SDS).

STORAGE AND STABILITY

Store at 2-8 °C.

Components are stable once opened until the expiry date marked in the label if they are stored well closed and care is taken to prevent contamination during their use.

On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 2 months.

Indications of deterioration: Absorbance of the blank below the limit indicated in "Test Parameters".

ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

□ Pyridoxal phosphate (BioSystems cod. 11666):

C. Reagent: Pyridoxal phosphate 10 mmol/L. 5 mL.

□ Biochemistry Calibrator (BioSystems cod. 18011) or Biochemistry Calibrator Human (BioSystems cod. 18044).

REAGENT PREPARATION

Reagents are provided ready to use.

PROCEDURE		
Reading mode	monochromatic	monochromatic
Main filter	340	340
Reference filter	-	-
Sample	16.5	16.5
Vol. R1	160	160
Vol. R2	40	40
Reading 1 (cycle)	23	46
Reading 2 (cycle)	33	66
Predilution factor	-	-

	distilled water experimental calibrator	distilled water experimental calibrator
Blank type		
Calibration mode		
Number of calibrators	1	1
Calibration curve	-	-
OPTIONS		
Blank absorbance limit	1.400	1.400
Kinetic blank limit	-	-
Linearity limit	500	500
Substrate depletion	0.5	0.5

CALIBRATION AND BLANK

M215311-05

BioSystems S.A. Costa Brava, 30.

08030 Barcelona (Spain)

03/2018

www.biosystems.es

Quality System certified according to
EN ISO 13485 and EN ISO 9001 standards

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- Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. Burtis CA, Ashwood ER, Bruns DE. WB Saunders Co, 2005.
- Sociedad Española de Química Clínica, Comité Científico, Comisión de Enzimas. Método recomendado para la determinación en rutina de la concentración catalítica de la alanina aminotransferasa en suero sanguíneo humano. *Quim Clin* 1987; 6: 241-244.
- IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 4. Reference procedure for the measurement of catalytic concentration of alanine aminotransferase. *Clin Chem Lab Med* 2002; 40: 718-724.
- IFCC reference procedures for measurement of catalytic concentrations of enzymes; corrigendum, notes and useful advice. *Clin Chem Lab Med* 2010; 48: 615-621.
- Gella FJ, Olivella T, Cruz Pastor M, Arenas J, Moreno R, Durban R and Gómez JA. A simple procedure for routine determination of aspartate aminotransferase and alanine aminotransferase with pyridoxal phosphate. *Clin Chim Acta* 1985; 153: 241-247.
- Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
- World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1, Rev.2; 2002.

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