12/2019

Barcelona, July, 20, 2020

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I, Laura López, Regulatory Affairs Manager at BioSystems S.A., hereby certify, that attached documents are true and accurate copies of User Manual of the product and translation of this document into Russian language in accordance with the documents.

Also, I hereby certify that the attached documents are for purposes of product registration in Russia and in the countries of Eurasian Economic Union.

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Laura López  
Regulatory Affairs Manager  
BIOSYSTEMS S.A.

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PROTEIN (TOTAL)

COD 21513 2 x 60 mL + 2 x 20 mL	COD 23513 1 x 60 mL + 1 x 20 mL
Only for <i>in vitro</i> use in the clinical laboratory	

PROTEIN (TOTAL)  
BIURET

INTENDED USE

Reagents for the measurement of protein (total) concentration in human serum or plasma. The obtained values are useful as an aid in the evaluation of the states of dehydration and water intoxication and salt retention syndromes, and to control the evolution of multiple myeloma, and to control Wadenström macroglobulinemia.

This reagent is for use in the BioSystems BA analyzers.

CLINICAL SIGNIFICANCE

Most of the plasma proteins are synthesized by the liver. The major exception to this is the immunoglobulins which are produced by plasma cells found in the spleen, lymph nodes and bone marrow.

The two general causes of alterations of serum total protein are a change in the volume of plasma water and a change in the concentration of one or more of the serum proteins.

Hyperproteinemia can be caused by dehydration (inadequate water intake, severe vomiting, diarrhea, Addison's disease, diabetic acidosis) or as a result of an increase in the concentration of specific proteins (immunoglobulins in chronic infections, multiple myeloma)<sup>1,2</sup>.

Hypoproteinemia may be caused by hemodilution (salt retention syndromes, massive intravenous infusions), by an impaired synthesis (severe malnutrition, chronic liver disease, intestinal malabsorptive disease), or by an excessive protein loss due to a chronic kidney disease or severe burns<sup>1,2</sup>.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE OF THE METHOD

Protein in the sample reacts with copper (II) ion in alkaline medium forming a coloured complex that can be measured by spectrophotometry<sup>3</sup>.

CONTENTS

	COD 21513	COD 23513
A. Reagent	2 x 60 mL	1 x 60 mL
B. Reagent	2 x 20 mL	1 x 20 mL

COMPOSITION

A. Reagent. Sodium hydroxide 0.4 mol/L, sodium tartrate 90 mmol/L.

**DANGER: H315: Causes skin irritation. H318: Causes serious eye damage. P280: Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P332+P313: If skin irritation occurs: Get medical advice/attention. P362: Take off contaminated clothing and wash before reuse.**

B. Reagent. Sodium hydroxide 0.4 mol/L, sodium tartrate 90 mmol/L, Copper (II) acetate 21 mmol/L, potassium iodide 60 mmol/L.

**DANGER: H315: Causes skin irritation. H318: Causes serious eye damage. P280: Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P332+P313: If skin irritation occurs: Get medical advice/attention. P362: Take off contaminated clothing and wash before reuse.**

STORAGE AND STABILITY

Store at 2-30°C.

Components are stable once opened until the expiry date marked in the label if they are stored well closed and care is taken to prevent contamination during their use.

On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 42 days.

Indications of deterioration: Absorbance of the blank over the limit indicated in "Test Parameters".

WARNING AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

Biochemistry Calibrator (BioSystems cod. 18011) or Biochemistry Calibrator Human (BioSystems cod. 18044).

REAGENT PREPARATION

Reagents are provided ready to use.

SAMPLES

Serum or heparinized plasma collected by standard procedures. Stable for 4 weeks at 4-8°C. Anticoagulants other than heparin should not be used.

CALIBRATION

A reagent blank should be done every day and a calibration at least every 21 days, after reagent lot change or as required by quality control procedures.

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, 18009 and 18042) and II (cod. 18007, 18010 and 18043) to verify the accuracy of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits.

REFERENCE VALUES

Serum, adults<sup>1</sup>:

Ambulatory	64-83 g/L
Recumbent	60-78 g/L

Concentrations are lower in child. Plasma total protein concentration is 2 to 4 g/L higher due to the presence of fibrinogen as well as some other trace proteins<sup>1</sup>.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

The metrological characteristics described below have been obtained using a BA400 analyzer and following the guidelines of the Clinical & Laboratory Standards Institute (CLSI).

- Detection limit: 0.800 g/L.
- Linearity limit: 150 g/L. For samples with higher values, dilute manually or refer to the Test Parameterization for Automatic dilution (note that all these samples will be diluted with the same dilution ratio).
- Precision:

Mean concentration	Repeatability (CV)	Within-laboratory (CV)
50.0 g/L	0.5 %	1.6 %
81.8 g/L	0.6 %	1.1 %

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

LIMITATIONS OF THE PROCEDURE

- Interferences: bilirubin (up to 30 mg/dL), hemolysis (hemoglobin up to 500 mg/dL) and lipemia (triglycerides up to 975 mg/dL) do not interfere. Other drugs and substances may interfere<sup>3</sup>.

BIBLIOGRAPHY

1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. Burtis CA, Ashwood ER, Bruns DE. WB Saunders Co, 2005.
2. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AAC Press, 2001.
3. Gornall AG, Bardwell CS, David MM. Determination of serum proteins by means of the Biuret reaction. *J Biol Chem* 1949; 177: 751-766.
4. World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations Document WHO/DIL/LAB/99.3, Rev 2, 2002. <http://www.who.int/iris/handle/10665/65957>.
5. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AAC Press, 2000.

TEST PARAMETERS

R1: use Reagent A.

R2: use Reagent B.

	BA200	BA400
<b>GENERAL</b>	<b>PROTEIN TOTAL</b>	<b>PROTEIN TOTAL</b>
Name	T-PROT	T-PROT
Short name	serum / plasma	serum / plasma
Sample type	differential bioreagent	differential bioreagent
Analysis mode		
Unit	g/L	g/L
Decimals	0	0
Reaction type	increasing	increasing
<b>PROCEDURE</b>		
Reading mode	monochromatic	monochromatic
Main filter	535	535
Reference filter	-	-
Sample	4	4
Vol. R1	225	225
Vol. R2	75	75
Reading 1 (cycle)	17	35
Reading 2 (cycle)	34	67
Predilution factor	-	-
Reduced factor	2	2
<b>CALIBRATION AND BLANK</b>		
Blank type	distilled water	distilled water
Calibration mode	experimental calibrator	experimental calibrator
Number of calibrators	1	1
Calibration curve	-	-
<b>OPTIONS</b>		
Blank absorbance limit	0,150	0,150
Kinetic blank limit	-	-
Linearity limit	150	150
Substrate depletion	-	-

**Content modifications**

1. Protein (Total) (BA) 2x60+2x60 ml
2. Protein (Total) (AU) 2x60+2x60 ml
3. Protein (Total) (D) 2x60+2x60 ml
4. Protein (Total) (F) 2x60+2x60 ml
5. Protein (Total) (BA) 1x60+1x60 ml
6. Protein (Total) (AU) 1x60+1x60 ml
7. Protein (Total) (D) 1x60+1x60 ml
8. Protein (Total) (F) 1x60+1x60 ml

**DISPOSAL AND DESTRUCTION OF MEDICAL DEVICES**

National standards must be followed when transferring and disposing of the product, its components and packaging.

Chemicals and finished and used reagents must be disposed of in accordance with regional regulations.

**LIABILITY AND GUARANTEES**

Any change to the set and/or exchange or replacement of any components of different batches from one set to another may negatively affect the expected results and the validity of the analysis as a whole. In the event of such changes and/or substitutions, any requirement to replace the set is invalid.

Claims submitted due to incorrect interpretation of the analysis results by the client are also invalid.

However, in case of any claims, the manufacturer's liability should not exceed the cost of the kit. Damage caused to the kit during transportation is not the responsibility of the manufacturer. The manufacturer's warranties of BioSystems S. A. apply to the fact that this product meets the specifications stated on the label, but only if the product is used by qualified personnel in accordance with the instructions set out on the label and the materials are used before the expiration date indicated on the label.

**MAINTENANCE**

No requirement

**MAINTENANCE SERVICE**

Not required

**POTENTIAL CONSUMER**

Professional user: a laboratory specialist who has the appropriate work permit

**INFORMATION ABOUT THE PRESENCE OF A MEDICINAL PRODUCT IN A MEDICAL DEVICE**

Does not contain

**INDICATIONS FOR USE**

Not applicable

**CONTRAINDICATIONS TO USE**

Not applicable

**STERILITY**

It is delivered in non-sterile form. Sterility is not required

**INFORMATION ABOUT THE POSSIBILITY OF REUSE**

Not re-usable

**SYMBOLS and LABELS**

Symbol	Name
	See manual
	Catalogue number
	Medical device <i>in vitro</i>
	Lot
	European Community symbol
	Temperature
	Expire before
	Manufacturer

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**AUTORISED REPRESENTATIVE of the MANUFACTURER**

“MD-Consulting and Development” LLC, 129337, Russian Federation, Moscow, Yaroslavskoye Shosse, 6/2, office 106

**MANUFACTURER AND ADDRESS OF PRODUCTION**

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1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. Burtis CA, Ashwood ER, Bruns DE. WB Saunders Co, 2005.
2. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
3. Gornall AG, Bardawill CS, David MM. Determination of serum proteins by means of the Biuret reaction. J Biol Chem 1949; 177: 751-766.
4. World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations Document WHO/DIL/LAB/99.1, Rev.2; 2002. <http://www.who.int/iris/handle/10665/65957>.
5. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

**BIOSYSTEMS, S.A.**  
Costa Brava, 30,  
08030 BARCELONA  
Tel. +34 933 11 00 00

Regulatory Affairs Manager  
BIOSYSTEMS S.A.  
Laura Lopez

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