12/2019

FE3575200

Barcelona, July, 20, 2020

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I, Laura López, Regulatory Affairs Manager at BioSystems S.A., hereby certify, that attached documents are true and accurate copies of User Manual of the product and translation of this document into Russian language in accordance with the documents.

Also, I hereby certify that the attached documents are for purposes of product registration in Russia and in the countries of Eurasian Economic Union.

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Laura López  
Regulatory Affairs Manager  
BIOSYSTEMS S.A.

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# ANTI-STREPTOLYSIN O (ASO)

COD 22923 2 x 60 mL + 2 x 15 mL	COD 23923 1 x 60 mL + 1 x 15 mL
Only for <i>in vitro</i> use in the clinical laboratory	

# ANTI-STREPTOLYSIN O (ASO) LATEX

### INTENDED USE

Reagents for the measurement of anti-streptolysin O (ASO) concentration in human serum. The obtained values are useful as an aid in the diagnosis and monitoring of streptococcal infections.

This reagent is for use in the BioSystems BA analyzers or in other analyzer with similar performance characteristics.

### CLINICAL SIGNIFICANCE

Anti-streptolysin O are the specific antibodies to streptolysin O, an extracellular enzyme produced by Lancefield group A,  $\beta$ -hemolytic streptococci (*Streptococcus pyogenes*). Antibodies against streptolysin O can be detected from one week to one month after the onset of a streptococcal infection. *Streptococcus pyogenes* causes a wide variety of upper respiratory infections such as acute pharyngitis. Other manifestations of *Streptococcus pyogenes* infection include glomerulonephritis, rheumatic fever, bacterial endocarditis and scarlet fever<sup>1,4</sup>.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

### PRINCIPLE OF THE METHOD

Serum anti-streptolysin O (ASO) causes agglutination of latex particles coated with streptolysin O. The agglutination of the latex particles is proportional to the ASO concentration and can be measured by turbidimetry<sup>5</sup>.

### CONTENTS

	COD 22923	COD 23923
A. Reagent	2 x 60 mL	1 x 60 mL
B. Reagent	2 x 15 mL	1 x 15 mL

### COMPOSITION

A. Reagent: Tris buffer 20 mmol/L, sodium chloride 150 mmol/L, sodium azide 0.95 g/L, pH 8.2.  
B. Reagent: Suspension of latex particles coated with streptolysin O, sodium azide 0.95 g/L (Note 1).

### STORAGE AND STABILITY

Store at 2-8°C.

Components are stable once opened until the expiry date marked in the label if they are stored well closed and care is taken to prevent contamination during their use.

On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 2 months.

Indications of deterioration: Absorbance of the blank over the limit indicated in "Test Parameters".

### ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

S. ASO Standard: 1 x 1 mL (BioSystems cod. 31119), Human serum. ASO concentration is given on the label. The concentration value is traceable to the Biological Reference Material WHO 94/572 (National Institute for Biological Standards and Control, NIBSC).

*Human serum used in the preparation of the standard has been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, the standard should be handled cautiously as potentially infectious.*

Reconstitute with 1.00 mL of distilled water. Stable for 30 days at 2-8°C.

### REAGENT PREPARATION

Reagents are provided ready to use.

### SAMPLES

Serum collected by standard procedures.

Anti-streptolysin O in serum is stable for 7 days at 2-8°C.

### CALIBRATION

A reagent blank should be done every day and a calibration at least every 2 months, after reagent lot change or as required by quality control procedures.

### QUALITY CONTROL

It is recommended to use the Rheumatoid Control Serum level I (cod. 31213) and II (cod. 31214) to verify the accuracy of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits.

### REFERENCE VALUES

Serum<sup>6</sup>

Adults: < 200 IU/mL

Children: < 150 IU/mL

This range is given for orientation only; each laboratory should establish its own reference range.

### METROLOGICAL CHARACTERISTICS

The metrological characteristics described below have been obtained using a BA400 analyzer and following the guidelines of the Clinical & Laboratory Standards Institute (CLSI).

- Detection limit: 8.4 IU/mL.

- Linearity limit: 800 IU/mL.

- Precision:

Mean concentration	Repeatability (CV)	Within-laboratory (CV)
187 IU/mL	1.8 %	3.2 %
255 IU/mL	1.8 %	3.0 %

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

### LIMITATIONS OF THE PROCEDURE

- Interferences: bilirubin (up to 20 mg/dL) and hemolysis (hemoglobin up to 1000 mg/dL), lipemia (triglycerides up to 1000 mg/dL) and rheumatoid factors (up to 2200 IU/mL) do not interfere. Other drugs and substances may interfere<sup>7</sup>.

- Zone effect: falsely low values are obtained when ASO is present in the sample at a concentration higher than 4000 IU/mL.

### NOTES

1. Shake the Reagent B vial gently before using.

### BIBLIOGRAPHY

1. Bisno AL. Group A streptococcal infections and acute rheumatic fever. *N Engl J Med* 1991; 325: 783-93.
2. Stevens DL. Invasive group A streptococcal disease. *Clin Infect Dis* 1992; 14: 2-11.
3. Immunology and Serology in Laboratory Medicine, 2<sup>nd</sup> edition. Turgeon mL, Mosby, 1996.
4. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
5. Borque L, Rus A, Dubois H. Automated determination of streptolysin O antibodies by turbidimetric latex immunoassay method. *J Clin Immunoassay* 1992; 15: 182-6.
6. Klein GC, Baker CN, Jones WL. Upper limits of normal antistreptolysin O and antideoxyribonuclease B titers. *Appl Microbiol* 1971; 21: 758-60.
7. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

### TEST PARAMETERS

These reagents may be used in several automatic analyzers. Specific instructions for application in many of them are available on request.

R1: use Reagent A, R2: use Reagent B.

	BA200	BA400
<b>GENERAL</b>		
Name	ASO	ASO
Short name	ASO	ASO
Sample type	serum	serum
Analysis mode	endpoint bireagent	endpoint bireagent
Unit	IU/mL	IU/mL
Decimals	0	0
Reaction type	increasing	increasing
<b>PROCEDURE</b>		
Reading mode	monochromatic	monochromatic
Main filter	535	535
Reference filter	-	-
Sample	3	3
Vol. R1	240	240
Vol. R2	60	60
Reading 1 (cycle)	34	67
Reading 2 (cycle)	-	-
Predilution factor	-	-
<b>CALIBRATION AND BLANK</b>		
Blank type	distilled water	distilled water
Calibration mode	experimental calibrator	experimental calibrator
Number of calibrators	1	1
Calibration curve	-	-
<b>OPTIONS</b>		
Blank absorbance limit	0.900	0.900
Kinetic blank limit	-	-
Linearity limit	800	800
Substrate depletion	-	-

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1. Bisno AL. Group A streptococcal infections and acute rheumatic fever. *Engl J Med* 1991; 325: 783-93.
2. Stevens DL. Invasive group A streptococcal disease. *Clin Infect Dis* 1992; 14: 2-11.
3. Immunology and Serology in Laboratory Medicine, 2<sup>nd</sup> edition. Turgeon mL. Mosby, 1996.
4. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
5. Borque L, Rus A, Dubois H. Automated determination of streptolysin O antibodies by turbidimetric latex immunoassay method. *J Clin Immunoassay* 1992; 15: 182-6.
6. Klein GC, Baker CN, Jones WL. Upper limits of normal antistreptolysin O and antideoxyribonuclease B titers. *Appl Microbiol* 1971; 21: 758-60.
7. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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Regulatory Affairs Department  
BIOSYSTEMS S.A.

**BIOSYSTEMS, S.A.**  
Costa Brava, 30,  
08030, BARCELONA  
Tel. +34 93 311 00 00

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